

ARG82225 Bile Acid Assay Kit (Fluorometric)

Package: 100 assay
Store at: -80°C

Summary


Product Description	ARG82225 Bile Acid Assay Kit (Fluorometric) is a detection kit for the quantification of Bile Acid in serum, plasma, cell lysate and tissue lysate.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Bile Acid
Conjugation Note	λ ex/em = 560 / 590 nm (fluorometric)
Sensitivity	1 μ M
Sample Type	Serum, plasma, cell lysate and tissue lysate.
Standard Range	0.39 - 25 μ M
Sample Volume	50 μ l

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	1 hour

Properties

Form	Liquid
Storage instruction	Store the kit at -80°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82225 Bile Acid Assay Kit (Fluorometric) standard curve image

ARG82225 Bile Acid Assay Kit (Fluorometric) results of a typical standard.