

ARG82281 Glucose Assay Kit

Package: 100 tests
Store at: 4°C, -20°C

Summary


Product Description	ARG82281 Glucose Assay Kit is a detection kit for the quantification of Glucose in serum, plasma, cell culture supernatants, food and beverage.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Glucose
Conjugation Note	Read at 630 nm.
Sensitivity	0.7 mg/dl (38.85 μ M)
Detection Range	0.7 - 300 mg/dl (38.85 μ M - 16.6 mM)
Sample Type	Serum, plasma, cell culture supernatants, food and beverage.
Sample Volume	As little as 5 μ l samples

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	~ 15 min

Properties

Form	Liquid
Storage instruction	Store components at 4°C or -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82281 Glucose Assay Kit standard curve image

ARG82281 Glucose Assay Kit results of a typical data with optical density reading at 630 nm.