

## ARG82304 Human HAV IgM antibody ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG82304 Human HAV IgM antibody ELISA Kit is an Enzyme Immunoassay kit for the determination of Human HAV IgM antibody in serum and plasma (EDTA, heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Human HAV IgM antibody
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	61.7 GBU/ml
Sample Type	Serum and plasma (EDTA, heparin, citrate).
Standard Range	Cut-off
Sample Volume	100 µl
Precision	Intra-Assay CV: < 15% Inter-Assay CV: < 20%

### Application Instructions

Assay Time	~ 3 hours
------------	-----------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.