

ARG82305 HBV surface antigen / HBsAg ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

| | |
|---------------------|--|
| Product Description | ARG82305 HBV surface antigen / HBsAg ELISA Kit is an Enzyme Immunoassay kit for the determination of HBsAg in serum and plasma (EDTA, heparin, citrate). |
| Tested Reactivity | HBV |
| Tested Application | ELISA |
| Target Name | HBV surface antigen / HBsAg |
| Conjugation | HRP |
| Conjugation Note | Substrate: TMB and read at 450 nm. |
| Sensitivity | 0.2 ng/ml |
| Sample Type | Serum and plasma (EDTA, heparin, citrate). |
| Standard Range | Cut-off |
| Sample Volume | 50 µl |
| Precision | Intra-Assay CV: 7.5% Inter-Assay CV: 25.7% |

Application Instructions

| | |
|------------|-----------|
| Assay Time | ~ 2 hours |
|------------|-----------|

Properties

| | |
|---------------------|--|
| Form | 96 well |
| Storage instruction | Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |