

ARG82556 Human GDNF Receptor alpha 3 ELISA Kit

Package: 96 wells
Store at: 4°C

Component

Cat. No.	Component Name	Package	Temp
ARG82556-001	Antibody-coated microplate	8 X 12 strips	4°C. Unused strips should be sealed tightly in the air-tight pouch.
ARG82556-002	Standard	2 X 10 ng/vial	4°C
ARG82556-003	Standard/Sample diluent	30 ml (Ready to use)	4°C
ARG82556-004	Antibody conjugate concentrate (100X)	1 vial (100 µl)	4°C
ARG82556-005	Antibody diluent buffer	12 ml (Ready to use)	4°C
ARG82556-006	HRP-Streptavidin concentrate (100X)	1 vial (100 µl)	4°C
ARG82556-007	HRP-Streptavidin diluent buffer	12 ml (Ready to use)	4°C
ARG82556-008	25X Wash buffer	20 ml	4°C
ARG82556-009	TMB substrate	10 ml (Ready to use)	4°C (Protect from light)
ARG82556-010	STOP solution	10 ml (Ready to use)	4°C
ARG82556-011	Plate sealer	4 strips	Room temperature

Summary

Product Description	ARG82556 Human GDNF Receptor alpha 3 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human GDNF Receptor alpha 3 in serum, plasma (EDTA, heparin, citrate) and cell culture supernatants.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	GDNF Receptor alpha 3
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	23.45 pg/ml
Sample Type	Serum, plasma (EDTA, heparin, citrate) and cell culture supernatants.
Standard Range	46.9 - 3000 pg/ml
Sample Volume	100 µl
Precision	Intra-Assay CV: 5.3%

Inter-Assay CV: 6.9%

Alternate Names

GDNFR-alpha-3; GFR-alpha-3; GDNF family receptor alpha-3; GDNFR3; GDNF receptor alpha-3

Application Instructions

Assay Time ~ 5 hours


Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	GFRA3
Gene Full Name	GDNF family receptor alpha 3
Background	The protein encoded by this gene is a glycosylphosphatidylinositol(GPI)-linked cell surface receptor and a member of the GDNF receptor family. It forms a signaling receptor complex with RET tyrosine kinase receptor and binds the ligand, artemin. [provided by RefSeq, Jul 2008]
Function	Receptor for the glial cell line-derived neurotrophic factor, ARTN (artemin). Mediates the artemin-induced autophosphorylation and activation of the RET receptor tyrosine kinase. [UniProt]
Highlight	Related products: GDNF Receptor antibodies; GDNF Receptor ELISA Kits; New ELISA data calculation tool: Simplify the ELISA analysis by GainData
PTM	N-glycosylated. [UniProt]
Cellular Localization	Cell membrane; Lipid-anchor, GPI-anchor. [UniProt]

Images



ARG82556 Human GDNF Receptor alpha 3 ELISA Kit standard curve image

ARG82556 Human GDNF Receptor alpha 3 ELISA Kit results of a typical standard run with optical density reading at 450 nm.