

## ARG82658 Formate Assay Kit (Colorimetric)

Package: 100 assay  
Store at: RT, -20°C, -80°C

### Summary


Product Description	ARG82658 Formate Assay Kit (Colorimetric) is a detection kit for the quantification of Formate in serum, plasma, urine, cell lysate and cell culture supernatants.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Formate
Conjugation	Un-conjugated
Conjugation Note	Read at 450 nm (colorimetric)
Sensitivity	12.5 µM
Sample Type	Serum, plasma, urine, cell lysate and cell culture supernatants.
Standard Range	12.5 - 800 µM
Sample Volume	50 µl

### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	~ 1 hour

### Properties

Form	Liquid
Storage instruction	Store components at RT, -20°C or -80°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82658 Formate Assay Kit (Colorimetric) standard curve image

ARG82658 Formate Assay Kit (Colorimetric) results of a typical standard run with optical density reading at 450 nm, using colorimetric reading.