

## ARG82660 Gentamicin ELISA Kit

Package: 96 wells  
Store at: 4°C, -20°C

### Summary

Product Description	ARG82660 Gentamicin ELISA Kit is an Enzyme Immunoassay kit for the quantification of Gentamicin in serum, plasma (heparin, citrate), cell lysate, tissue homogenates and food samples.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	Gentamicin
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	23 nM
Sample Type	Serum, plasma (heparin, citrate), cell lysate, tissue homogenates and food samples.
Standard Range	23 - 1500 nM
Sample Volume	50 µl

### Application Instructions


Assay Time	~ 3 hours
------------	-----------

### Properties

Form	96 well
Storage instruction	Store components at 4°C or -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Highlight	Related products:  New ELISA data calculation tool: <a href="#">Simplify the ELISA analysis by GainData</a>
-----------	--



ARG82660 Gentamicin ELISA Kit standard curve image

ARG82660 Gentamicin ELISA Kit results of a typical standard run with optical density reading at 450 nm.