

## ARG82791 Human 2-MET (Metanephrine + Normetanephrine) ELISA Kit


Package: 2x96 wells

Store at: 4°C

### Summary


Product Description	ARG82791 Human 2-MET (Metanephrine + Normetanephrine) ELISA Kit is an Enzyme Immunoassay kit for the quantification of 2-MET (Metanephrine + Normetanephrine) in Human plasma (EDTA, heparin).
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	

### Images



ARG82791 Human 2-MET (Metanephrine + Normetanephrine) ELISA Kit standard curve image

ARG82791 Human 2-MET (Metanephrine + Normetanephrine) ELISA Kit results of a typical standard run with optical density reading at 450 nm.



ARG82791 Human 2-MET (Metanephrine + Normetanephrine) ELISA Kit standard curve image

ARG82791 Human 2-MET (Metanephrine + Normetanephrine) ELISA Kit results of a typical standard run with optical density reading at 450 nm.