


## ARG82793 Human Normetanephrine ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG82793 Human Normetanephrine ELISA Kit is an Enzyme Immunoassay kit for the quantification of Normetanephrine in Human urine.
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	

### Images



ARG82793 Human Normetanephrine ELISA Kit standard curve image

ARG82793 Human Normetanephrine ELISA Kit results of a typical standard run with optical density reading at 450 nm.