


ARG82796
Human 17-OH Progesterone ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG82796 Human 17-OH Progesterone ELISA Kit is an Enzyme Immunoassay kit for the quantification of 17-OH Progesterone in Human serum and plasma (EDTA).
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	

Images



ARG82796 Human 17-OH Progesterone ELISA Kit standard curve image

ARG82796 Human 17-OH Progesterone ELISA Kit results of a typical standard run with optical density reading at 450 nm.