

## ARG82799 Human Cortisol ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG82799 Human Cortisol ELISA Kit provides competitive immunoenzymatic colorimetric method for the quantification of free Cortisol concentration in Human urine.
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	Analytical specificity: No interference has been found for Albumin, Acetylsalicylic Acid, Ibuprofen and Ascorbic Acid.
Target Name	Cortisol
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	0.47 ng/ml
Sample Type	Urine
Standard Range	1 - 200 ng/ml
Sample Volume	10 µl
Precision	Intra-Assay CV: 7.4% Inter-Assay CV: 10.8%

### Application Instructions

Assay Time	1 hour 15 min
------------	---------------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Highlight	Related products: <a href="#">Cortisol antibodies</a> ; <a href="#">Cortisol ELISA Kits</a> ; New ELISA data calculation tool: <a href="#">Simplify the ELISA analysis by GainData</a>
Research Area	Metabolism kit; Signaling Transduction kit