

ARG82805 Human DHEA (sulfate form) ELISA Kit

Package: 96 wells

Store at: 4°C

Summary

Product Description	ARG82805 Human DHEA (sulfate form) ELISA Kit is an Enzyme Immunoassay kit for the quantification of DHEA (sulfate form) in Human saliva.
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	Cross-Reactivity: DHEA: 100% DHEA-S: 90% Androsterone-S-Na: 48% Androstendione: 20% Etiocolanone-S-Na: 0.2% 5-Androstendione, Testosterone, Progesterone, 17-OH Progesterone and Estrone: 0.01% Cortisol and Cholesterol: 0.001%
Target Name	DHEA (sulfate form)
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	0.05 ng/ml
Sample Type	Saliva
Standard Range	0.2 - 12 ng/ml
Sample Volume	50 µl
Precision	Intra-Assay CV: < 7.8% Inter-Assay CV: < 14.9%

Application Instructions

Assay Time	30 min
------------	--------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.