

## ARG82836 Lactate (D-Lactate) Assay Kit

Package: 100 tests  
Store at: -20°C

### Summary


Product Description	ARG82836 Lactate (D-Lactate) Assay Kit is a detection kit for the quantification of Lactate (D-Lactate) in serum, plasma and cell culture media.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Lactate (D-Lactate)
Conjugation Note	Read at 565 nm.
Sensitivity	0.05 mM
Detection Range	0.2 - 2 mM
Sample Type	Serum, plasma and cell culture media.
Sample Volume	20 µl
Alternate Names	Lactic acid; D-lactic acid

### Application Instructions

Assay Time	~ 20 min
------------	----------

### Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82836 Lactate (D-Lactate) Assay Kit standard curve image

ARG82836 Lactate (D-Lactate) Assay Kit results of a typical standard run with optical density reading at 565 nm.