

ARG82851 Mouse FABP2 / Intestinal FABP ELISA Kit

Package: 96 wells
Store at: 4°C

Component

Cat. No.	Component Name	Package	Temp
ARG82851-001	Antibody-coated microplate	8 X 12 strips	4°C. Unused strips should be sealed tightly in the air-tight pouch.
ARG82851-002	Standard	2 X 10 ng/vial	4°C
ARG82851-003	Standard/Sample diluent	30 ml (Ready to use)	4°C
ARG82851-004	Antibody conjugate concentrate (100X)	1 vial (100 µl)	4°C
ARG82851-005	Antibody diluent buffer	12 ml (Ready to use)	4°C
ARG82851-006	HRP-Streptavidin concentrate (100X)	1 vial (100 µl)	4°C
ARG82851-007	HRP-Streptavidin diluent buffer	12 ml (Ready to use)	4°C
ARG82851-008	25X Wash buffer	20 ml	4°C
ARG82851-009	TMB substrate	10 ml (Ready to use)	4°C (Protect from light)
ARG82851-010	STOP solution	10 ml (Ready to use)	4°C
ARG82851-011	Plate sealer	4 strips	Room temperature

Summary

Product Description	ARG82851 Mouse FABP2 / Intestinal FABP ELISA Kit is an Enzyme Immunoassay kit for the quantification of Mouse FABP2 / Intestinal FABP in serum, plasma (EDTA, heparin) and cell culture supernatants.
Tested Reactivity	Ms
Tested Application	ELISA
Target Name	FABP2 / Intestinal FABP
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	15.6 pg/ml
Sample Type	Serum, plasma (EDTA, heparin) and cell culture supernatants.
Standard Range	31.2 - 2000 pg/ml
Sample Volume	100 µl
Precision	Intra-Assay CV: 5.9%

Inter-Assay CV: 7.0%

Alternate Names

Intestinal-type fatty acid-binding protein; I-FABP; Fatty acid-binding protein, intestinal; Fatty acid-binding protein 2; FABPI

Application Instructions

Assay Time ~ 5 hours

Properties

Form 96 well

Storage instruction Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.

Note For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol FABP2


Gene Full Name fatty acid binding protein 2, intestinal

Background The protein encoded by this gene is an intracellular fatty acid-binding protein that participates in the uptake, intracellular metabolism, and transport of long-chain fatty acids. The encoded protein is also involved in the modulation of cell growth and proliferation. This protein binds saturated long-chain fatty acids with high affinity, and may act as a lipid sensor to maintain energy homeostasis. [provided by RefSeq, Aug 2017]

Function FABP are thought to play a role in the intracellular transport of long-chain fatty acids and their acyl-CoA esters. FABP2 is probably involved in triglyceride-rich lipoprotein synthesis. Binds saturated long-chain fatty acids with a high affinity, but binds with a lower affinity to unsaturated long-chain fatty acids. FABP2 may also help maintain energy homeostasis by functioning as a lipid sensor. [UniProt]

Cellular Localization Cytoplasm. [UniProt]

Images



ARG82851 Mouse FABP2 / Intestinal FABP ELISA Kit standard curve image

ARG82851 Mouse FABP2 / Intestinal FABP ELISA Kit results of a typical standard run with optical density reading at 450 nm.