

## ARG82865 Human Aspergillus fumigatus IgM Antibody ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

|                     |  |
|---------------------|--|
| Product Description | ARG82865 Human Aspergillus fumigatus IgM Antibody ELISA Kit has been designed for the qualitative determination of specific IgM antibodies against Aspergillus fumigatus in serum and plasma (Citrate, heparin). |
| Tested Reactivity   | Hu   |
| Tested Application  | ELISA  |
| Target Name         | Aspergillus fumigatus IgM antibody   |
| Sensitivity         | Diagnostic Sensitivity: 91 %   |
| Sample Type         | Serum and plasma (Citrate, heparin).   |
| Standard Range      | Cut-off  |
| Sample Volume       | 5-10 µl  |

### Application Instructions

|            |                 |
|------------|-----------------|
| Assay Time | 1 h, 30, 15 min |
|------------|-----------------|

### Properties

|                     |  |
|---------------------|--|
| Form                | 96 well  |
| Storage instruction | Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note                | For laboratory research only, not for drug, diagnostic or other use.   |

### Bioinformation

|               |   |
|---------------|---|
| Research Area | Microbiology and Infectious Disease kit |
|---------------|---|