

## ARG82885 Human Measles IgG antibody ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

|                     |  |
|---------------------|--|
| Product Description | ARG82885 Human Measles IgG antibody ELISA Kit has been designed for the qualitative determination of specific IgG antibodies against Measles in serum and plasma (Citrate, heparin). |
| Tested Reactivity   | Hu   |
| Tested Application  | ELISA  |
| Target Name         | Measles IgG antibody   |
| Sensitivity         | Diagnostic Sensitivity: 97 %   |
| Sample Type         | Serum and plasma (Citrate, heparin).   |
| Standard Range      | Cut-off  |
| Sample Volume       | 5-10 µl  |

### Application Instructions

|            |                 |
|------------|-----------------|
| Assay Time | 1 h, 30, 15 min |
|------------|-----------------|

### Properties

|                     |  |
|---------------------|--|
| Form                | 96 well  |
| Storage instruction | Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note                | For laboratory research only, not for drug, diagnostic or other use.   |

### Bioinformation

|               |  |
|---------------|--|
| Research Area | Immune System kit; Microbiology and Infectious Disease kit |
|---------------|--|