

## ARG82949 Monkey IgM ELISA kit

Package: 96 wells  
Store at: 4°C

### Summary

|                     |   |
|---------------------|---|
| Product Description | ARG82949 Monkey IgM ELISA kit is an Enzyme Immunoassay kit for the quantification of Monkey IgM in serum, plasma and cell culture supernatants. |
| Tested Reactivity   | Mk  |
| Tested Application  | ELISA   |
| Target Name         | IgM   |
| Conjugation         | HRP   |
| Conjugation Note    | Substrate: TMB and read at 450 nm.  |
| Sensitivity         | 1.56 ng/ml  |
| Sample Type         | Serum, plasma and cell culture supernatants.  |
| Standard Range      | 3.125 - 200 ng/ml   |
| Sample Volume       | 100 µl  |

### Application Instructions

|            |             |
|------------|-------------|
| Assay Time | ~ 3.5 hours |
|------------|-------------|

### Properties

|                     |  |
|---------------------|--|
| Form                | 96 well  |
| Storage instruction | Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note                | For laboratory research only, not for drug, diagnostic or other use.   |