

## ARG82956 Mouse IgG2b ELISA kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG82956 Mouse IgG2b ELISA kit is an Enzyme Immunoassay kit for the quantification of Mouse IgG2b in serum, plasma and cell culture supernatants.
Tested Reactivity	Ms
Tested Application	ELISA
Target Name	IgG2b
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	0.156 ng/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	0.313 - 20 ng/ml
Sample Volume	100 µl

### Application Instructions

Assay Time	~ 3.5 hours
------------	-------------

### Properties

Form	96 well
Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.