

ARG82974 Human BMP4 ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG82974 Human BMP4 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human BMP4 in serum, plasma, cell culture supernatant samples.
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	Cross-Reactivity: Not reacts with Human BMPR-IB/FChimera, BMP5, BMP6, BMP7 and BMP8 proteins.
Target Name	BMP4
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	7.8 pg/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	15.6 - 1000 pg/ml
Sample Volume	100 µl
Alternate Names	BMP2B; BMP-4; OFC11; Bone morphogenetic protein 2B; BMP2B1; Bone morphogenetic protein 4; MCOPS6; ZYME; BMP-2B

Application Instructions

Assay Time ~ 4 hours

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	BMP4
Gene Full Name	bone morphogenetic protein 4
Background	The protein encoded by this gene is a member of the bone morphogenetic protein family which is part of the transforming growth factor-beta superfamily. The superfamily includes large families of growth and differentiation factors. Bone morphogenetic proteins were originally identified by an ability of demineralized bone extract to induce endochondral osteogenesis in vivo in an extraskelatal site. This particular family member plays an important role in the onset of endochondral bone formation in humans, and a reduction in expression has been associated with a variety of bone diseases, including the heritable disorder Fibrodysplasia Ossificans Progressiva. Alternative splicing in the 5' untranslated

region of this gene has been described and three variants are described, all encoding an identical protein. [provided by RefSeq, Jul 2008]

Function

Induces cartilage and bone formation. Also act in mesoderm induction, tooth development, limb formation and fracture repair. Acts in concert with PTHLH/PTHRP to stimulate ductal outgrowth during embryonic mammary development and to inhibit hair follicle induction (By similarity). [UniProt]

Highlight

Related products:

[BMP antibodies](#); [BMP ELISA Kits](#); [BMP recombinant proteins](#);

New ELISA data calculation tool:

[Simplify the ELISA analysis by GainData](#)