

ARG83059 Toxoplasma gondii IgA ELISA

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG83059 Toxoplasma gondii IgA ELISA is an Enzyme Immunoassay kit for the qualitative determination of Toxoplasma gondii IgA in human serum or plasma (citrate, heparin).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Toxoplasma gondii IgA
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	diagn. 91.7 %
Detection Range	Cut-off
Sample Type	serum, plasma (citrate, heparin)
Sample Volume	100 µl

Application Instructions

Assay Time	~2 hour
------------	---------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Full Name	Toxoplasma gondii
Background	Toxoplasma gondii is an obligate intracellular parasitic protozoan (specifically an apicomplexan) that causes toxoplasmosis. Found worldwide, T. gondii is capable of infecting virtually all warm-blooded animals, but felids are the only known definitive hosts in which the parasite may undergo sexual reproduction.
Highlight	