

## ARG83068 Echinococcus / hydatid disease antibody ELISA Kit

Package: 96 wells

Store at: 4°C

### Summary

Product Description	ARG83068 Echinococcus / hydatid disease antibody ELISA Kit is an enzyme immunoassay kit for the quantification Echinococcus antibody in animals serum.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	Echinococcus
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Detection Range	Cut-off
Sample Type	serum
Standard Range	Cut-off
Sample Volume	100 µl
Alternate Names	Echinococcus, hydatid disease

### Application Instructions

Assay Time	~2 hour
------------	---------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Highlight