

ARG83099 Benzonase Nuclease ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	Benzonase Nuclease ELISA Kit is designed for the quantitative detection of residual nuclease content in intermediates, semifinished products and finished products of various biological products by using a double-antibody sandwich method.
Tested Application	ELISA
Sensitivity	0.05 ng/mL
Standard Range	0.05-5 ng/mL
Sample Volume	100 µL


Application Instructions

Assay Time	~3 hours
------------	----------

Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Images



ARG83099 Benzonase Nuclease ELISA Kit standard curve image

ARG83099 Benzonase Nuclease ELISA Kit results of standard run with optical density reading at 450 nm