

ARG83104 Mycoplasma DNA Detection Kit

Package: 100 tests
Store at: -20°C

Summary

Product Description	Mycoplasma DNA Detection Kit is designed for the detection of mycoplasma contamination in master cell bank, working cell bank, cells for clinical use and biological products. Mycoplasma DNA Detection Kit conforms to relevant regulations about mycoplasma testing in EP2.6.7 and JP XVI. Mycoplasma DNA Detection Kit adopts the qPCR-fluorescent probe method. The kit is a rapid, specific and reliable device and can finish the detection within 2 hours.
Tested Application	FuncSt

Application Instructions

Assay Time	~ 2 hour
------------	----------


Properties

Storage instruction	Store the kit at -20°C at all times. All reagents must be kept on ice during the entire experiment.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Highlight

Images



ARG83104 Mycoplasma DNA Detection Kit Standard curve (Positive)

ARG83104 Mycoplasma DNA Detection Kit results of standard run after 48 cycle



ARG83104 Mycoplasma DNA Detection Kit Standard curve (Negative)

ARG83104 Mycoplasma DNA Detection Kit results of standard run after 48 cycle