


ARG83115 Monkey Activin A ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG83115 Monkey Activin A ELISA Kit is an Enzyme Immunoassay kit for the quantification of Monkey Activin A in Serum, Plasma, Saliva, Cell culture supernatants.
Tested Reactivity	Mk
Tested Application	ELISA
Specificity	There is no detectable cross-reactivity with other relevant proteins.
Target Name	Activin A
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	

Images



ARG83115 Monkey Activin A ELISA Kit standard curve image

ARG83115 Monkey Activin A ELISA Kit ARG83115 results of a typical standard run with optical density reading at 450 nm.