


ARG83129 Porcine BMP2 ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG83129 Porcine BMP2 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Porcine BMP2 in Serum, Bone tissue, Cell culture supernatants.
Tested Reactivity	Pig
Tested Application	ELISA
Specificity	There is no detectable cross-reactivity with other relevant proteins.
Target Name	BMP2
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	

Images



ARG83129 Porcine BMP2 ELISA Kit standard curve image

ARG83129 Porcine BMP2 ELISA Kit results of a typical standard run with optical density reading at 450 nm.