


## ARG83130 Monkey BMP2 ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

|                     |   |
|---------------------|---|
| Product Description | ARG83130 Monkey BMP2 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Monkey BMP2 in Serum, Bone tissue, Cell culture supernatants. |
| Tested Reactivity   | Mk  |
| Tested Application  | ELISA   |
| Specificity         | There is no detectable cross-reactivity with other relevant proteins.   |
| Target Name         | BMP2  |
| Conjugation         | HRP   |
| Conjugation Note    | Substrate: TMB and read at 450 nm.  |
| Sensitivity         |   |

### Images



ARG83130 Monkey BMP2 ELISA Kit standard curve image

ARG83130 Monkey BMP2 ELISA Kit results of a typical standard run with optical density reading at 450 nm.