


## ARG83131 Rabbit BMP2 ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG83131 Rabbit BMP2 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Rabbit BMP2 in Serum, Bone tissue, Cell culture supernatants.
Tested Reactivity	Rb
Tested Application	ELISA
Specificity	There is no detectable cross-reactivity with other relevant proteins.
Target Name	BMP2
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	

### Images



ARG83131 Rabbit BMP2 ELISA Kit standard curve image

ARG83131 Rabbit BMP2 ELISA Kit results of a typical standard run with optical density reading at 450 nm.