


## ARG83188 Monkey IGFBP7 ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG83188 Monkey IGFBP7 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Monkey IGFBP7 in Serum, Plasma, Urine and Cell culture supernatants.
Tested Reactivity	Mk
Tested Application	ELISA
Specificity	There is no detectable cross-reactivity with other relevant proteins.
Target Name	IGFBP7
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	

### Images



ARG83188 Monkey IGFBP7 ELISA Kit standard curve image

ARG83188 Monkey IGFBP7 ELISA Kit results of a typical standard run with optical density reading at 450 nm.