

## ARG83363 Listeria monocytogenes p60 ELISA Kit

Package: 96 wells  
Store at: 4°C, -20°C

### Summary

Product Description	ARG83363 Listeria monocytogenes p60 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Listeria monocytogenes p60 in Food, Fecal and Environmental
Tested Reactivity	Bacteria
Tested Application	ELISA
Specificity	This ELISA kit is design for Listeria monocytogenes p60 protein
Target Name	Listeria monocytogenes
Conjugation	HRP
Conjugation Note	Read at 450 nm
Sensitivity	25 ng/ml
Sample Type	Food, Fecal and Environmental
Standard Range	0.625 - 40 ng/mL
Sample Volume	100 µL

### Application Instructions


Assay Time	~5 hours
------------	----------

### Properties

Form	96 well
Storage instruction	Store components at 4°C or -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Microbiology and Infectious Disease antibody
---------------	--



ARG83363 Listeria monocytogenes p60 ELISA Kit standard curve image

ARG83363 Listeria monocytogenes p60 ELISA Kit results of a typical standard run with optical density reading at 450 nm.

---