

ARG83379 Cas9 ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C, -80°C

Summary


Product Description	ARG83379 Cas9 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Cas9 in Cell Lysate and Tissue Lysate
Tested Reactivity	Other
Tested Application	ELISA
Target Name	Cas9
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	1 ng/mL
Sample Type	Cell Lysate and Tissue Lysate
Standard Range	1.56 - 100 ng/mL
Sample Volume	100 µL

Application Instructions

Assay Time	~3 hours
------------	----------

Properties

Form	96 well
Storage instruction	Store components at 4°C, -20°C and -80°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG83379 Cas9 ELISA Kit standard curve image

ARG83379 Cas9 ELISA Kit results of a typical standard run with optical density reading at 450 nm.