

ARG83386 Adenovirus Hexon ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG83386 Adenovirus Hexon ELISA Kit is an Enzyme Immunoassay kit for the quantification of Adenovirus Hexon in Various biological samples types.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	Adenovirus Hexon
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	0.1 ng/mL
Sample Type	Various biological samples types
Standard Range	0.156 - 10 ng/mL
Sample Volume	100 µL

Application Instructions

Assay Time	~5 hours
------------	----------

Properties

Form	96 well
Storage instruction	Store components at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.