

ARG83389
Human PCDH17 / PCDH68 ELISA KitPackage: 96 wells
Store at: 4°C**Summary**

Product Description	ARG83389 Human PCDH17 / PCDH68 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human PDGF AB in serum, plasma (heparin, EDTA) and cell culture supernatants.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	PCDH17 / PCDH68
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	