

ARG83524 Integral Reagent Kit (ELISA Development Kit)

Package: 1 kit(5 plates)

Component

| Cat. No. | Component Name | Package | Temp |
|-------------|--------------------------------------|---------------------------|--------------------------|
| ARG83524-01 | 96 well microplates | 40 X 12 strips (5 plates) | Room temperature |
| ARG83524-02 | PBS | 60mL | 4°C |
| ARG83524-03 | Blocking buffer | 2X 65mL | 4°C |
| ARG83524-04 | 10X Standard / Sample Diluent Buffer | 20mL | 4°C |
| ARG83524-05 | Antibody Diluent Buffer | 3X 60mL | 4°C |
| ARG83524-06 | 10X Wash Buffer | 3X 60mL | 4°C |
| ARG83524-07 | Substrate Solution | 5X 12mL | 4°C (protect from light) |
| ARG83524-08 | Stop Solution | 1X 60mL | 4°C |
| ARG83524-09 | Plate Sealers | 10 adhesive strips | Room temperature |

Summary

Product Description