

**ARG83554**  
**Lysozyme Activity Assay Kit**Package: 100 tests  
Store at: 4°C

### Summary

Product Description	ARG83554 Lysozyme Activity Assay Kit can be used to measure Lysozyme in Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Lysozyme
Conjugation Note	Read at 450 nm
Sensitivity	50 U/ml
Sample Type	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Standard Range	100 U/ml - 20000 U/ml
Sample Volume	10 µl
Alternate Names	muramidase, N-acetylmuramide glycanhydrolase

### Application Instructions

Assay Time	5 min
------------	-------

### Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.