

ARG83560 Butyrylcholinesterase Activity Assay Kit

Package: 100 tests

Store at: 4°C

Summary

Product Description	ARG83560 Butyrylcholinesterase Activity Assay Kit can be used to measure Butyrylcholinesterase in Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Butyrylcholinesterase
Conjugation Note	Read at 405 nm
Sensitivity	0.05 µmol/ml
Sample Type	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Standard Range	0.1 µmol/ml - 5 µmol/ml
Sample Volume	10 µl
Alternate Names	EC 3.1.1.8; BChE; Acylcholine acylhydrolase

Application Instructions

Assay Time	2 min
------------	-------

Properties

Storage instruction Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.

Note For laboratory research only, not for drug, diagnostic or other use.