

ARG83627 Human IL22 ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG83627 Human IL22 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human IL22 in serum, plasma and cell culture supernatants.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	IL22
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	3.90625 pg/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	7.8125 - 500 pg/ml
Sample Volume	50 µl
Precision	