

ARG83631
Mouse IL22 ELISA KitPackage: 96 wells
Store at: 4°C, -20°C**Summary**

Product Description	ARG83631 Mouse IL22 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Mouse IL22 in serum, plasma and cell culture supernatants.
Tested Reactivity	Ms
Tested Application	ELISA
Target Name	IL22
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	3.90625 pg/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	7.8125 - 500 pg/ml
Sample Volume	50 µl
Precision	